

VITAPM[®]9

Beauty. Strength. Versatility.



VITAPM[®]9 is VITA's latest restorative material in the VM line of small particle, 3D-Master shaded porcelain systems with excellent wear rates. PM9 is an extremely versatile press-to-zirconia system that enables technicians to consistently produce beautiful, high-strength restorations.

Pressable VITAPM[®]9 is based on the proven fine-structure ceramic VITAVM[®]9 and is used for pressing to yttria-stabilized ZrO₂ frameworks in the CTE range of 10.4 - 10.6 x 10⁻⁶, such as VITA In-Ceram YZ and other zirconium oxide substructure materials for natural looking high-strength anterior and posterior crowns and bridges.

VITAPM[®]9 may also be pressed without the support of a zirconia substructure for the fabrication of single or multi-surface inlays, onlays, partial crowns, veneers and anterior crowns using the staining/glazing or layering technique. Restorations made from VITAPM[®]9 can be easily characterized using the layering technique with the entire range of VITAVM[®]9 ceramic materials

or individualized (staining technique) with VITA Akzent or VITA Shading Paste stains. Each of the ten shaded pellets in the system is available in two levels of translucency.

Shades Include:

0M1P-O 0M2P-O 1M1P-O 1M2P-O 2M1P-O
2M2P-O 2M3P-O 3M1P-O 3M2P-O 3M3P-O



0M1P-T 0M2P-T 1M1P-T 1M2P-T 2M1P-T
2M2P-T 2M3P-T 3M1P-T 3M2P-T 3M3P-T



VITA[®]



800-828-3839

800-263-4778 in Canada • www.vident.com

VIDENT[™]

The Pellets

VITAPM[®]9 two gram pellets are available in ten VITA 3D-Master shades, each with two different levels of translucency.

The shade of the Opaque pellet of VITAVM[®]9 has been matched with YZ Coloring Liquid. The opaque version is primarily suitable for the press-to technique and the translucent version is suitable for the staining and layering technique for unsupported pressing.



Available in packs containing 5 pellets in 10 shades: OM1P, OM2P, 1M1P, 1M2P, 2M1P, 2M2P, 2M3P, 3M1P, 3M2P, 3M3P each available in O (opaque) and T (translucent) variations.

VITAPM[®]9 Complete Assortment

- 20 packs (5pcs. each) of VITAPM[®]9 press pellets (2g) in 10 shades and 2 translucency levels (opaque (O) and translucent (T) in VITA 3D-Master shades.
- 1 pack of 50 disposable press plungers
- **PM9** Investment (56x100g)
- **PM9** Investment Liquid (900ml)
- 1 Investment System (ring, formers), 3 pieces (200g)
- 1 syringe (12ml) of VITA Firing Paste
- 1 **PM9** Shade Guide (Opacues shades)

| Part No. | Product Description: |
|-----------|--|
| VX36-9001 | PM9 Kit |
| VX36-9201 | PM9 Pellets, 0M1 P-T, 2g/ea Translucent, 5pk |
| VX36-9202 | PM9 Pellets, 0M2 P-T, 2g/ea Translucent, 5pk |
| VX36-9211 | PM9 Pellets, 1M1 P-T, 2g/ea Translucent, 5pk |
| VX36-9212 | PM9 Pellets, 1M2 P-T, 2g/ea Translucent, 5pk |
| VX36-9221 | PM9 Pellets, 2M1 P-T, 2g/ea Translucent, 5pk |
| VX36-9222 | PM9 Pellets, 2M2 P-T, 2g/ea Translucent, 5pk |
| VX36-9223 | PM9 Pellets, 2M3 P-T, 2g/ea Translucent, 5pk |
| VX36-9231 | PM9 Pellets, 3M1 P-T, 2g/ea Translucent, 5pk |
| VX36-9232 | PM9 Pellets, 3M2 P-T, 2g/ea Translucent, 5pk |
| VX36-9233 | PM9 Pellets, 3M3 P-T, 2g/ea Translucent, 5pk |
| VX36-9301 | PM9 Pellets, 0M1 P-O, 2g/ea Opaque, 5pk |
| VX36-9302 | PM9 Pellets, 0M2 P-O, 2g/ea Opaque, 5pk |
| VX36-9311 | PM9 Pellets, 1M1 P-O, 2g/ea Opaque, 5pk |
| VX36-9312 | PM9 Pellets, 1M2 P-O, 2g/ea Opaque, 5pk |
| VX36-9321 | PM9 Pellets, 2M1 P-O, 2g/ea Opaque, 5pk |
| VX36-9322 | PM9 Pellets, 2M2 P-O, 2g/ea Opaque, 5pk |
| VX36-9323 | PM9 Pellets, 2M3 P-O, 2g/ea Opaque, 5pk |
| VX36-9331 | PM9 Pellets, 3M1 P-O, 2g/ea Opaque, 5pk |
| VX36-9332 | PM9 Pellets, 3M2 P-O, 2g/ea Opaque, 5pk |
| VX36-9333 | PM9 Pellets, 3M3 P-O, 2g/ea Opaque, 5pk |
| J36-001 | PM9 Investment System, 3 parts |
| J36-002 | PM9 Investment, 56 x 100g |
| J36-003 | PM9 Investment Liquid, 900ml |
| J36-004 | PM9 Disposable Plungers (50) |
| A36-001 | PM9 Firing Paste, 12g |
| A36-002 | PM9 Firing Paste, 3x12g |
| GX36-901 | PM9 Color Indicator |

VITA[®]

VIDENT[™]

800-828-3839

800-263-4778 in Canada • www.vident.com